

## ACCESSORIES EVE IN/TR/NEO

Please submit the request-form via Email ([sales@stephan-gmbh.com](mailto:sales@stephan-gmbh.com)) or Fax (06439/9125-111) to our sales department. You will receive your offer as soon as possible.

### SINGLE-USE



Description	Art.-No.	Amount
Breathing circuit, adult, emergency, Length 180 cm	1 070 61 120	
Length 240 cm	1 070 61 140	
Length 300 cm	1 070 61 141	
Length 440 cm	1 070 61 142	



Description	Art.-No.	Amount
Breathing circuit, adult, intensive	1 070 61 125	



Description	Art.-No.	Amount
Breathing circuit, pediatrics, with integrated flow sensor	1 070 61 124	

## ACCESSORIES EVE IN/TR/NEO

### SINGLE-USE



Description	Art.-No.	Amount
Breathing circuit neonatal (distal) incl. expiration valve , Length 180 cm	1 070 61 189	



Description	Art.-No.	Amount
Flow sensor, adult, dead space 11 ml, flow range $\geq 20$ l/min, Length 200 cm	1 070 61 127	
Length 260 cm	1 070 61 144	
Length 320 cm	1 070 61 145	
Length 460 cm	1 070 61 146	



Description	Art.-No.	Amount
Flow sensor, pediatric, dead space 2.7 ml, flow range $< 20$ l/min	1 070 61 128	



Description	Art.-No.	Amount
Respiration sensor (Graseby Capsule) for external trigger control, EVE NEO, Length 150 cm	1 035 60 103	
Length 200 cm	1 035 60 203	

## ACCESSORIES EVE IN/TR/NEO


### SINGLE-USE



Description	Art.-No.	Amount
Expiration valve (distal)	1 070 61 072	



Description	Art.-No.	Amount
Adapter pressure measurement NCPAP	1 070 61 137	



Description	Art.-No.	Amount
Duplex hose with PNT-adaptor, Length 200 cm	1 070 61 138	

## ACCESSORIES EVE IN/TR/NEO

### REUSABLE



Description	Art.-No.	Amount
Breathing circuit emergency for adults, autoclavable	1 070 61 190	



Description	Art.-No.	Amount
PNT-B, flow sensor neonatal, 0 to +/- 12 l/min, Length 200 cm	1 070 61 129	



Description	Art.-No.	Amount
Breathing circuit EVE NEO, autoclavable	1 070 61 197	

Company

Address

Person in charge

Email

Phone

Date

Signature