

ACCESSORIES EVE IN/TR/NEO, EVA, EVA NEO

Please submit the request-form via Email (sales@stephan-gmbh.com) or Fax (+49 6439/9125-118) to our sales department. You will receive your offer as soon as possible.

SINGLE-USE



Description	Art.-No.	Amount
Breathing circuit, adult, emergency, Length 180 cm	1 070 61 120	
Length 240 cm	1 070 61 140	
Length 300 cm	1 070 61 141	
Length 440 cm	1 070 61 142	



Description	Art.-No.	Amount
Breathing circuit, adult, intensive	1 070 61 125	



Description	Art.-No.	Amount
Breathing circuit for EVE/EVA with expiration valve distal, pressure measuring connection and flow sensor	1 070 61 198	
Breathing circuit for EVE with expiration valve distal, pressure measuring connection and flow sensor	1 070 61 162	



Description	Art.-No.	Amount
Breathing circuit, pediatrics, with integrated flow sensor	1 070 61 124	

ACCESSORIES EVE IN/TR/NEO, EVA, EVA NEO

SINGLE-USE




Description	Art.-No.	Amount
Breathing circuit neonatal (distal) incl. expiration valve , Length 180 cm	1 070 61 189	



Description	Art.-No.	Amount
Flow sensor, adult, dead space 11 ml, flow range ≥ 20 l/min, Length 200 cm	1 070 61 127	
Length 260 cm	1 070 61 144	
Length 320 cm	1 070 61 145	
Length 460 cm	1 070 61 146	



Description	Art.-No.	Amount
Flow sensor, pediatric, dead space 2.7 ml, flow range < 20 l/min	1 070 61 128	



Description	Art.-No.	Amount
Flow-Sensor, Neonatology, dead space 1.3 ml, Length 200 cm	1 070 61 158	

ACCESSORIES EVE IN/TR/NEO, EVA, EVA NEO

SINGLE-USE



Description	Art.-No.	Amount
Expiration valva distal with pressure measuring connection	1 070 61 081	




Description	Art.-No.	Amount
Expiration valve "Combi", Adult Flow sensor incl.	1 070 61 072	



Description	Art.-No.	Amount
Expiration valve "Combi", Paediatrics Flow sensor incl.	1 070 61 077	

ACCESSORIES EVE IN/TR/NEO, EVA, EVA NEO

SINGLE-USE



Description	Art.-No.	Amount
Duplex hose with PNT-adaptor, Length 200 cm	1 070 61 138	



Description	Art.-No.	Amount
Respiration sensor (Graseby Capsule) for external trigger control, EVE NEO, Length 150 cm	1 035 60 103	
Length 200 cm	1 035 60 203	



Description	Art.-No.	Amount
① Flow sensor electronic, Neonatology	1 070 61 150	
② Flow sensor electronic, Adult	1 070 61 152	

ACCESSORIES EVE IN/TR/NEO, EVA, EVA NEO

REUSABLE



Description	Art.-No.	Amount
PNT-B, flow sensor neonatology, 0 to +/- 12 l/min, Length 200 cm	1 070 61 129	



Description	Art.-No.	Amount
PNT-D, flow sensor Adult, 0 to +/- 80 l/min, Length 200 cm	1 070 61 147	



Description	Art.-No.	Amount
① Flow sensor electronic, Neonatology	1 760 61 151	
② Flow sensor electronic, Adult	1 760 61 153	



Description	Art.-No.	Amount
Exchangeable battery EVE/EVA	1 070 61 037	

Company _____ Address _____ Person in charge _____

Email _____ Phone _____

Date _____ Signature _____